

MD300 C3 Pulse Oximeter

- Displays SpO₂, Pulse Rate and SpO₂ wavelength
- Multi-directional screen display mode
- Audible Pulse rate
- Displays signal strength
- Auto-off function
- 12 month replacement warranty



About the MD300 C3

The MD300 C3 is a great addition to the line of pulse oximeters offered by Maxtec. Its compact design allows it to be used for pediatric and adult applications in many different areas such as hospital, homecare and EMS. This innovative technology is an excellent source for the inexpensive, accurate pulse oximeter that you have been looking for.

Specifications

Patient Range:	Adult, Pediatric and Neonatal
SpO ₂ :	Range: 0%-99%, Resolution: 1%, Accuracy: 80%-99%: ±2%, 70%-80%: ±3%, 0%-69%: unspecified
Heart (Pulse) Rate:	Range: 0-254 bpm, Resolution: 1bpm, Accuracy: ±2bpm or ±2%, Voice: PR tone modulation, Information: Signal strength, bargraph
Alarm:	Alarm Finger out, Low power indicator
Operating Temperature:	5°C - 40°C
Storage Temperature:	-10°C - 40°C
Operating Humidity:	15%RH - 80%RH
Storage Humidity:	10%RH - 80%RH
Mechanicals:	Size: 58 x 32 x 28mm, Weight: 50g (including 2 batteries)
Power:	Type: 2AAA alkaline batteries, Operation time: about 30 hours for normal operation
Display:	Type: LCD, Parameters: SpO ₂ , PR, SpO ₂ waveform, Mode: 6 display modes, Brightness: adjustable 10 levels

Part Numbers

MD300 C3 Pulse OximeterR204P19