

• PULSOX-2

- **Lightweight (2 ounces), easy to use**
- **80 hrs of use on 2 AAA batteries**
- **Spot check**
- **Low cost**
- **2 year warranty with optional Maxtec® protection plan**



Specifications

Measuring method:	Dual-wavelength pulse-type oximeter
Measuring range:	SpO2: 0 to 100% Pulse rate: 20 to 250 bpm
Accuracy:	SpO2: ± 2% (70 to 100%, 1 s.d.) Pulse rate: ± 2 bpm
Display:	Oxygen saturation (SpO2) Pulse rate number Pulse level meter Error messages
Equipment Classification:	Protection against electric shock: Internally powered Type of applied Part: BF Splash-proof equipment: not designed to be submerged into water. Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide. Mode of operation of Equipment: Continuous while in Use (IEC60601-1)
Operation temperature/humidity range:	0 to 40°C (32 to 104°F) 30 to 95% relative humidity with no condensation
Storage temperature/humidity range:	-10 to 60°C (14 to 140°F) 30 to 95% relative humidity with no condensation
Power:	2 AAA-size batteries
Dimensions (W x H x D):	69 x 60 x 28 mm (2-11/16 x 2-3/8 x 1-1/8 in.)
Weight (without batteries):	57 g (2 oz.)
Accessories:	Standard Accessory; Neck strap

Part Numbers

Pulsox-2 R204P12