

Coiled Hose

- Easy to move with equipment
- DISS Fittings (handtight and hextight)
- Self coiling prevents touching the ground
- Meets NFPA 99 Standard for manufactured assemblies



15' Single Coil Oxygen Hose



15' Single Coil Air Hose



Cross Section

A Solution for Gas Delivery

Maxtec® offers unique oxygen and air self-coiling hoses. They are made of a high-grade polyurethane material to help keep the work space clean and tidy. They are available in both single and dual coils with DISS fittings only.

Specifications

Burst Rate:	1,000 psi
Connectors:	DISS Male, Female
Green Hose:	Oxygen
Yellow Hose:	Air

Part Numbers

15' Single Coil Oxygen Hose (DISS fittings)	R127P35	10' Dual Coil Blender Hose (DISS fittings)	R129P01
15' Single Coil Air Hose (DISS fittings).	R128P34	15' Dual Coil Ventilator Hose (DISS fittings).	R129P02
50' Extension Tubing	R132P01		

Great for home patient use (Image not shown)