

MaxVenturi™

- Provides users with the ability to mix ambient air & oxygen w/o the need of piped medical grade air
- Long battery life on analyzer (approximately 5,000 hours)
- 24 month warranty
- Independent adjustable O2% and flow
- All-in-one oxygen/air mixer, MAXO2®+ analyzer, flow meter & built-in pole clamp



Patent Pending

The MaxVenturi™ is specially designed to interface with the high-flow humidification systems. It provides users with the ability to mix ambient air and oxygen without the need of piped medical grade air and creates the perfect system for patients needing high flow therapy. This unit also includes the Maxtec® analyzer technology as well as a high-flow flow meter, giving you all three features you desire in a high flow delivery system.

Part Numbers

MaxVenturi™ (with Max-250E Sensor)	R211P03-002
MaxVenturi™ (with Max-250E & NIST fittings)	R211P03-015