

MD300 C5 Pediatric Oximeter

- Specially designed to fit a wide range of pediatric patients
- Displays pulse rate and oxygen saturation
- Powered by 2 AAA batteries
- One button on/off
- 12 month warranty



About the MD300 C5

The MD300C5 is a great addition to the line of pulse oximeters offered by Maxtec. Its compact design allows it to be used on a wide range of pediatric applications in many different areas such as hospital, homecare and EMS. This innovative technology is an excellent source for the inexpensive, accurate pulse oximeter that you have been looking for.

Specifications

Patient Range:	Pediatric and Infant patients
SpO ₂ :	Range: 70%-100%, Resolution: 1%, Accuracy: 80%-100%: ±2%, 70%-79%: ±3%, 0%-69%: unspecified
Heart (Pulse) Rate:	Range: 30-235 bpm, Resolution: 1bpm, Accuracy: 30 -100 ±2bpm, 101 - 235 ±2%
Alarm:	Battery-low indicator
Operating Temperature:	5°C - 40°C
Storage Temperature:	-20°C - 55°C
Operating Humidity:	RH ≤ 85%
Storage Humidity:	RH ≤ 85%
Dimensions:	Size: 50 x 28 x 28mm, Weight: 25g (without batteries)
Power:	Type: 2AAA alkaline or rechargeable batteries, Operation time: about 30 hours for normal operation
Display:	Type: Dual color OLED display, Parameters: SpO ₂ , PR, Pulse bar, Mode: 6 display modes, Brightness: adjustable 10 levels

Part Numbers

MD300 C5 Pediatric Pulse Oximeter R204P20