

•• Pulsox-300

- **ALL NEW** Soft Tip Probe
- **4000** spot measurements on 1 AAA battery
- **Reinforced design** for increased durability
- **Large, easy-to-read LED display**
- **User-friendly design** with automatic power-off function
- **12 month warranty**



Spot Check Solution

A new addition to our pulse oximetry line-up, the Pulsox-300 offers durability and comfort in this all new wrist model oximeter. This unit is best used during stress tests, physician to patient testing and for simple spot checking situations.

Specifications

Model Dimensions:	2.68" (W) x 2.28" (H) x 0.59" (D)
Model Weight:	56 grams (2 ounces)
Measuring Range:	SpO ₂ : 0 - 100%, Pulse rate: 30 to 230 bpm
Accuracy:	SpO ₂ : ±2% (70 - 100%), Pulse Rate: ±2 bpm (30 to 230 range)
Warning Functions:	SpO ₂ : Low SpO ₂ warning, Battery indication, Error message indication
Operating Conditions:	0 to 40°C (32 to 104°F); 30 to 85% relative humidity with no condensation 700 to 1060hPa (altitude: -400 to 3000m)
Storage Conditions:	-10 to 60°C (14 to 140°F); 10 to 95% relative humidity with no condensation 700 to 1060hPa (altitude: -400 to 3000m)
Power Requirements:	1 AAA Alkaline battery (included)
Included Accessories:	Wrist strap, finger probe, 1 AAA battery

Part Numbers

Pulsox-300 (no probe) R124P06	Pulsox-300 (w/finger probe) R204P16
Pulsox-300 (w/spot-check probe) R204P16-001	Pulsox-300 (w/soft tip probe) R204P16-002