


# SpiroQuant A+ Flow Sensor

- Fast response time
- Reliable
- Accurate low pressure drop
- Accurate temperature compensation
- DEHP free
- 100% tested and inspected



## SpiroQuant A+ Flow Sensor

The SpiroQuant A+ is for measuring volume flows in the field of anaesthetics/ventilation.

### Specifications

Measurement Principle	Constant temperature hot wire - anemometry
Application	Dräger [843735 Sensor]
Accuracy	± 8% reading
Added Resistance	< 2.5mbar
Pressure Range	± 100mbar
Connections	Dräger specific
Operation temperature	+15 to +40°C
Storage Temperature	-20 to +40°C

### Part Numbers

SpiroQuant A+ Flow Sensor. . . . . R122P01

† Maxtec is not affiliated with nor endorsed by Dräger Medical GmbH. Sensor compatibility has not been approved by Dräger.