

• The SmartHook IV Pole Adapter

- **Secure Multiple Poles**
- **Standard 6" Length** (Center to Center)
- **Custom Lengths Available**
- **Keeps Poles Together, Neat and Orderly**
- **Transport Multiple Poles With Ease**
- **Snaps Securely To 1" Diameter Poles**



Pictured to right: SmartHook application on SmartStack Premium Poles

About the SmartHook:

Maxtec is pleased to announce our all-new **SmartHook** IV Pole adapter. Do you ever have the need to transport more than one IV pole? It can be a cumbersome task, especially when these poles are loaded with heavy equipment. With our SmartStack IV Pole combined with the new SmartHook adapter, moving multiple IV Poles is done quickly and with ease! The standard SmartHook was specially designed to snap securely around 1" diameter poles and provide 6" of separation when connected. Do you have a need for a longer SmartHook for your application? No problem! We can fully customize the length to your specifications. Contact Maxtec Customer Service for more information.

Specifications:

Length:	6" Center to center of hook (Custom sizes available)
Fits Pole Size:	Fits 1" (25mm) diameter poles
Material:	Delrin

Part Numbers:

SmartHook R105P16