

FireSafe™ Nozzles

- Slows the rate of a fire developing
- Prevents fire from reaching upstream equipment
- Increases chance of the fire extinguishing itself
- Minimize fire damage and cost of restoration
- Increased patient safety
- Sold in pack of ten



FireSafe™ with DISS Fittings



Inline FireSafe™

About the FireSafe™ Nozzle:

FireSafe™ Nozzles fit directly to your existing oxygen flowmeter outlet and can limit the impact of an oxygen fire by automatically arresting the oxygen flow. Acting like a fuse, it automatically isolates the oxygen supply providing increased patient safety while minimizing the potential for fire disaster.

A Comparison of Possible Risk Controls

CONTROL	RESIDUAL RISK COMPONENT	AUTOMATIC OPERATION	POTENTIAL BENEFIT IN OXYGEN THERAPY FIRE		
			PREVENTS IGNITION	ABATES OR EXTINGUISHES	AIDS EVACUATION
User Training	Occurrence		Yes		
Smoke Alarm	Severity	Yes			Yes
Fire Extinguisher	Severity			Yes	Yes
Sprinkler System	Severity	Yes		Yes	Yes
FireSafe™	Severity	Yes		Yes	Yes

Part Numbers

FireSafe™ Nozzle **with DISS Fittings** R110P02
 FireSafe™ Nozzle **Inline with 2" Tubing** R400P16-001

FireSafe™ Nozzle **Inline** R400P16