


SmartStack Premium I.V. Stand

- Unique stacking design with taller outer tube
- Large, smooth rolling casters with 2 wheel locks
- Includes 4 fixed Ram-horn hooks
- Easy, one-handed height adjustment
- 48 month warranty



The pole height can be adjusted quickly & with ease



The SmartStack I.V. Stand from Maxtec just got even better!

The SmartStack Premium combines our specially designed pole base with the existing features of the basic model and more. The Premium model features a new pole with an increased outer tube height. Furthermore, this premium pole features a one-handed height adjustment knob that releases and locks into place quickly and with ease.

Specifications

Collapsed Height:	Pole 54", Pole with Base 61" (measured from floor to top of hook assembly)
Extended Height:	Pole Extended 86", Pole with Base 93" (measured from floor to top of hook assembly)
Casters:	Neoprene Tread, Chrome Steel, 3" QTY: 5, 2 Locking, 3 Non-Locking
Weight Capacity:	100lbs Total, 50lbs on Top Portion, 11lb Maximum One Hook
Wheel Base Diameter:	22"
Nesting Base:	4" (center to center)
Weight:	14lbs
Shipping Weight:	Base: 83lbs, Poles: 28lbs (six units)
Pole Diameter:	25mm (bottom portion) 18mm (upper portion)
Base Material:	Powder coated aluminum
Pole Material:	Stainless steel
Ram-Horn Hooks:	4 fixed
Restrictions:	Not MRI Compatible

Part Numbers

SmartStack Premium Single-Pack (1pk) R100P49-001
 SmartStack Premium Fourteen-Pack (14pk) . . . R100P49-014
 Basket R105P05
 Blickman Patient Handle R105P06-001
 I.V. Pole Organizer R105P15

SmartStack Premium Two-Pack (2pk) R100P49-002
 Single Tank Holder R105P04
 Dual Tank Holder R105P04-001
 Patient Handle R105P06